

Mobius FoamFighter™

External anti foam system for Mobius 3L bioreactor

Date: 10.2.2017

Version: 1.3

External AF System


This autonomous anti Foam System detects Foam in Bioreactors and applies anti foam media from a syringe. Because of its external arrangement and optical sensor, the system can also be added when the experiment is already running. The foam detection level and detection distance can be adjusted.

Features:

- Status LED
- External control and feedback signals for Applikon Easy Control
- Programmable dosing intervals
- Optical foam detection
- Adjustable foam tolerance level
- Adjustable detection distance

This version fits:

- Mobius 3L Single use Bioreactor
- 10 ml Syringe (BD)



ReseaTech GmbH | Pestalozzistrasse 16 | 3400 Burgdorf | Switzerland
www.reseatech.ch | info@reseatech.ch

This information corresponds to the current state of knowledge. Reseatech reserves the right to make technical changes. Liability for consequential damage resulting from the use of Reseatech products is excluded

 **reseatech**
more than dosing!